

Plasma

Components	Contents	Barcode #	Form #
FRESH FROZEN PLASMA	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18211	2800
FRESH FROZEN PLASMA / AUTOMATED PHERESIS	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18211	2804
FRESH FROZEN PLASMA / PHERESIS	___mL containing approx. ___mL of CP2D Anticoagulant. Store at -18 C. or colder.	18211	2827
FRESH FROZEN PLASMA / DOUBLE UNIT	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18211 CUSTOM	2840
FRESH FROZEN PLASMA / DOUBLE UNIT	___mL containing approx. ___mL of Trisodium Citrate. Store at -18 C. or colder. Collected by Apheresis Technology.	18211 CUSTOM	2841
FRESH FROZEN PLASMA / PHERESIS	___mL containing approx. ___mL of Sodium Citrate. Store at -18 C. or colder.	18211 CUSTOM	2828
FRESH FROZEN PLASMA / PHERESIS	___mL containing approx. ___mL of CPDA-1. Store at -18 C. or colder.	18211 CUSTOM	2838
FRESH FROZEN PLASMA / AUTOMATED PHERESIS	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18212	2823
FRESH FROZEN PLASMA	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18212 CUSTOM	2858
FRESH FROZEN PLASMA / PHERESIS	250 mL containing approx. 25 mL of 4% Sodium Citrate. Store at -18 C. or colder.	18212 CUSTOM	2825
FRESH FROZEN PLASMA / PHERESIS	250 mL containing approx. 25 mL of CP2D Anticoagulant. Store at -18 C. or colder.	18212 CUSTOM	2860
FRESH FROZEN PLASMA / AUTOMATED PHERESIS	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18213 CUSTOM	2857
FRESH FROZEN PLASMA / THAWED	___mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.	18231	9034
FRESH FROZEN PLASMA / THAWED	___mL from 500 mL CPD Whole Blood. Store at 1 to 6 C.	18231 CUSTOM	9046
FRESH FROZEN PLASMA	___mL containing approx. ___mL of ACD-A Anticoagulant. Store at -18 C. or colder.	18251 CUSTOM	2864
CRYO POOR FROZEN PLASMA	___mL from 450 mL CPD Whole Blood. Store at -18 C. or colder.	18401	9040
CRYO POOR FROZEN PLASMA	___mL from 450 mL CP2D Whole Blood. Store at -18 C. or colder	18401	9039

LEGEND: * Custom Accommodation: 1 roll or more required with minimum order of \$50.00

Plasma

Components	Contents	Barcode #	Form #
CRYO POOR FROZEN PLASMA	___mL from 500 mL CPD Whole Blood. Store at -18 C. or colder.	18401	9042
PLASMA	___mL from 450 mL CPD Whole Blood. Store at -18 C. or colder.	18401	9010
PLASMA	___mL from 450 mL CP2D Whole Blood. Store at -18 C. or colder.	18401	9011
PLASMA	___mL from 500 mL CPD Whole Blood. Store at -18 C. or colder.	18401	9012
PLASMA	___mL from 450 mL CPDA-1 Whole Blood. Store at -18 C. or colder.	18401	9015
PLASMA	___mL from 450 mL of CPDA-1, CPD, CP2D Whole Blood. Store at -18 C. or colder.	18401	9009
PLASMA / AHF REMOVED	___mL from 450 mL CPDA-1 Whole Blood. Store below -18 C. or colder.	18401	9018
PLASMA	___mL from 500 mL CP2D Whole Blood. Store at 1 to 6 C.	18401 CUSTOM	9013
PLASMA / AHF REMOVED	___mL from 450 mL CPD Whole Blood. Store at 1 to 6 C.	18401 CUSTOM	9017
PLASMA / AHF REMOVED	Approx. 225 mL from CPDA-1 Whole Blood. Store below -18 C.	18402	9020
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 450 mL CP2D Whole Blood. Store at -18 C. or colder.	18435	9023
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 500 mL CPD Whole Blood. Store at -18 C. or colder.	18435	9024
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 450 mL CPDA-1 Whole Blood. Store below -18 C. or colder.	18435	9028
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 450 mL CPD Whole Blood. Store at -18 C. or colder.	18435	9022
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 500 mL CPDA-1 Whole Blood. Store at -18 C. or colder.	18435 CUSTOM	9029
PLASMA / CRYOPRECIPITATE REDUCED	___mL from 500 mL CP2D Whole Blood. Store at -18 C. or colder.	18435 CUSTOM	9014
THAWED PLASMA	___mL from 500 mL CPD Whole Blood. Store at 1 to 6 C.	18499 CUSTOM	7107

LEGEND: * Custom Accommodation: 1 roll or more required with minimum order of \$50.00

Plasma

Components	Contents	Barcode #	Form #
THAWED PLASMA	___ mL from 500 mL CPD Whole Blood. Store at 1 to 6 C.	18500 CUSTOM	7108
THAWED PLASMA	___ mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.	18500 CUSTOM	7105
THAWED PLASMA	225 +/- 25 mL from 450 mL CPD Whole Blood. Store at 1 to 6 C.	18500 CUSTOM	7102
CRYO POOR FROZEN PLASMA / THAWED	___ mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.	18501	9037
LIQUID PLASMA	___ mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.	18501	9032
LIQUID PLASMA	___ mL from 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.	18501	9030
LIQUID PLASMA	___ mL from 450 mL CPD Whole Blood. Store at 1 to 6 C.	18501	9025
LIQUID PLASMA	___ mL from 500 mL CPD Whole Blood. Store at 1 to 6 C.	18501	9026
THAWED PLASMA	___ mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.	18501 CUSTOM	7103
THAWED PLASMA	___ mL from 500 mL ___ Whole Blood. Store at 1-6 C.	18501 CUSTOM	7111
THAWED PLASMA	___ mL from 500 mL CPD Whole Blood. Store at 1 to 6 C.	18502 CUSTOM	7104
THAWED FRESH / FROZEN PLASMA	___ mL from 450 mL CPD ___ Whole Blood. Store at 1 to 6 C.	18900	9035
THAWED FRESH / FROZEN PLASMA	___ mL from 500 mL CPD ___ Whole Blood. Store at 1 to 6 C.	18900 CUSTOM	9047
FRESH FROZEN PLASMA / THAWED, POOLED	Approx. ___ mL from ___ donors. Collected in CPD. Store at 1 to 6 C.	19411 CUSTOM	9051
FRESH FROZEN PLASMA / DONOR RETESTED / DIVIDED	___ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C. or colder.	47201 CUSTOM	7052
FRESH FROZEN PLASMA / DONOR RETESTED / DIVIDED	___ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C. or colder.	47202 CUSTOM	7053
FRESH FROZEN PLASMA / DONOR RETESTED / DIVIDED	___ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C. or colder.	47203 CUSTOM	7054

**HUMAN
LIQUID
PLASMA**

PLASMA

// or more required with

8. 6

___ mL from
450 mL CP2D
Whole Blood.
Store at -18° C. or colder.



FORM #9011

FORM #M-88

LEG
Rev